

# TIRF Labs Technical Note 3

## Compare TIRF Microscopy Geometries

### Stray Light Errors in Objective-type TIRF

This Technical Note analyzes errors caused by the stray light in TIRF microscopy and their interference with the TIRF effect. It compares the advantages and limitations of three popular TIRF microscopy geometries, discusses the sources of stray light in each of the geometries, and suggests ways for stray light mitigation. This Note will help you in choosing a TIRF geometry well-suited for your studies.


Total Internal Reflection Fluorescence (TIRF) has established itself as an indispensable tool for single molecule detection, cell membrane studies, and other areas of lifesciences. TIRF selectively excites only a ~100-nm layer of the specimen next to the surface, while the bulk of specimen is not illuminated and does not fluoresce. This superior spatial selectivity (in comparison, the confocal scheme excites ~1,000 nm) allows for the minimization of the background signal and providing the contrast sufficient for the detection of single molecules.

The phenomenon of total internal reflection occurs at the interface between two optical media with different refractive indices, e.g. glass/water. If the angle of incidence is larger than critical, the incident light reflects back into the glass and generates a profile of exponentially decaying intensity at the interface, termed the Evanescent Wave (EW) (Fig. 1). The critical angle is given by eq. (1) (Fig.1). For the glass/water interface the angle is 62 degrees. EW intensity is maximal at the surface and decays exponentially with distance - eq. (3) (Fig.1). The depth of penetration depends on the wavelength, the angle of incidence, and the ratio of refractive indices eq. (2). For a typical TIRF experiment with a glass/water interface, the depth of penetration is ~100 nm or less.

**TIRF Theory and Practice.** In theory, the intensity of EW exponentially decays with distance, as shown in Fig. 1. In practice, however, scatter, reflections, and refractions produce undesirable rays of light, collectively termed “stray light.” Stray light contaminates the exponential decay of EW, excites the bulk of specimen, and deteriorates the TIRF effect, as shown in Fig. 2.

**Sources of Stray Light.** All optical materials, to a certain extent, scatter light, and all surfaces and interfaces between optical parts reflect, refract, and scatter light. Due to the combination of these factors, the undesirable stray light is present in all practical systems. The intensity of stray light and its interference with the TIRF effect critically depend on the optical scheme used to implement TIRF. It also depends on the quality of the optical elements used in the system, including the quality of the materials themselves, the quality of manufacturing, especially surface quality, and the accumulation of dust particles and other external contaminations.

**TIRF Geometries.** There are three popular TIRF microscopy geometries: through-objective, prism-, and lightguide-based optical schemes (o-TIRF, p-TIRF, and Ig-TIRF, respectively). For a number of reasons, the reasons rather of business nature than the technical performance, Nikon, Olympus, Zeiss, and Leica have aggressively marketed o-TIRF geometry. Objective-TIRF employs expensive high NA objectives and other sophisticated optics. A typical o-TIRF microscopy system costs ~\$80,000 or more. On the other hand, TIRF geometries that demonstrate superior signal-to-background ratios, namely p-TIRF and Ig-TIRF, that are cost-efficient and flexible, were not offered commercially, until recently. In the past, each research group built p-TIRF and Ig-TIRF systems on their own, until 2010, when TIRF Labs started to market p-TIRF and Ig-TIRF systems. Since 2010, our customers generated unique TIRF data and demonstrated superior advantages of p- and Ig-TIRF geometries on a number of applications, including single molecule detection [1, 2]. The sections below describe potential sources of stray light in each of the geometries and outline ways to mitigate the stray light interferences.





Compare TIRF Microscopy Geometries  
Stray Light Errors in Objective-type TIRF


**Objective-type TIRF Geometry.** Fig. 3 shows the schematics of o-TIRF system. The main peculiar feature of o-TIRF, which might appear elegant at first glance, is the use of the emission path for delivering the excitation light to the glass/water interface. This feature has been inherited from epi-fluorescence. The o-TIRF scheme uses angles of incidence greater than the critical for glass/water interface > 63 degrees. Microscope objectives with a Numerical Aperture (NA) smaller than 1.38 do not support such angles. Therefore, the o-TIRF scheme depends on specialized high-NA objectives to deliver the excitation light at angles greater than critical. Fig. 3 illustrates numerous potential sources of the stray light in the o-TIRF scheme.

**Sources of stray light in o-TIRF.** Not surprisingly, significant interferences of stray light have been reported for o-TIRF geometry [3-8]. The stray light excites fluorophores in the bulk of the specimen and deteriorates the TIRF effect. In many instances of o-TIRF, the intensity of stray light is too large to calculate the depth of penetration. Two major sources of stray light were identified as (i) originating from the TIRF objective, and (ii) originating from the rest of the microscope optics [4, 5]. Only minor contributions were detected due to the scatter at the TIRF glass/water interface and at refractive-index boundaries within the specimen, including live cells [4, 5].

The first group of sources (i) is related to undesirable scatter and reflections inside the TIRF objective. The intense excitation light travels through multiple lenses and interfaces of the objective on its way to and from the TIRF surface, as shown in Fig. 3. The quality of the glass itself, which should minimally scatter light, and the surface quality of optical parts manufactured from the glass, are critically important for minimizing the intensity of stray light. To our knowledge, the systematic comparison of TIRF objectives from the standpoint of the intensity of stray light has not been performed yet. Our own tests on small number of TIRF objectives, the analysis of the literature, and reports of our customers and colleagues indicate that Nikon TIRF objectives demonstrate the smallest intensity of stray light. Perhaps, the better quality of Nikon lenses is related to the fact that only Nikon makes its own high-quality glass, while other firms rely on glass acquired from external sources.

The second group of stray light (ii) originates from the optics inside the microscope. A significant amount of stray light is generated at the dichroic mirror. Even a high-quality dichroic mirror scatters and transmits certain portion of light, which, in a perfect world, should be ideally reflected and blocked. The scatter of the excitation light at the dichroic mirror, as well as leaking of it to and through the emission filter, results in increased background.

In o-TIRF, the intensity of stray light changes unpredictably with the angle of incidence and XYZ coordinates. It increases with the amount of imperfections located on the path of the excitation light. Certain types of imperfections are distributed randomly, while other types exhibit more systematic patterns of their occurrence. If the angle of incidence increases, the depth of penetration and the intensity of EW excitation decreases, while the intensity of stray light remains the same or increases. If you are performing variable angle TIRF experiments using o-TIRF geometry, the effect of stray light should be carefully taken into account. In certain cases, the intensity of stray light is comparable with that of EW. In such cases, the depth of penetration calculated using eq. (3) (Fig.1) does not describe the intensity profile anymore and can mislead the interpretation of biological TIRF images [3-8].



Compare TIRF Geometries <-> Stray Light Causes Errors in Objective-type TIRF





## Compare TIRF Microscopy Geometries Stray Light Errors in Objective-type TIRF

**Dichroic mirror and emission filter.** The dichroic mirror is the central element of o-TIRF geometry; its quality is critically important for the objective-type TIRF. The scatter at the dichroic mirror appears to be the main source of the stray light. Additionally, the beam of excitation light reflected from the dichroic beamsplitter must be focused at the back focal plane of the objective. A significant focal shift or a change in the focal spot size caused by a bend of the dichroic mirror can make it difficult to achieve TIRF, especially if the microscope has a limited ability to adjust the collimation of the excitation beam. Credit should go to Chroma and Semrock engineers, who recently made great progress in improving technical performance of dichroic and bandpass filters. The companies increased the thickness of dichroic mirrors to keep their flatness, which is necessary for precision focusing in o-TIRF, minimized surface roughness, as well as the density of pinholes [9, 10].


**o-TIRF stray light mitigation.** First, it appears to be rational to explore the opportunity of using alternative TIRF geometries for your study. If your study dictates the use of o-TIRF scheme, select a TIRF objective with minimal rate of stray light. Use the best quality dichroic mirror and the emission filter. If your excitation light is not monochromatic, use an additional excitation filter to block undesirable lines in the excitation light. Finally, make sure that the dichroic mirror, emission filter and other accessible optical parts are free from dust particles and other contaminations.


**Prism-based TIRF.** p-TIRF geometry has been shown to provide the best signal-to-background ratio [1-3]. Indeed, p-TIRF is a scheme, where excitation and emission channels are naturally independent, as is illustrated in Fig. 4. The excitation light travels through the prism, reflects at the glass/water interface, and escapes through the opposite facet of the prism. Potential sources of the stray light are just a few reflections at the prism facets, prism/slide interface, and scatter inside the prism itself. These few sources of stray light are located away from the EW region and the detection channel. In p-TIRF, stray light minimally interferes with the EW region. Only a small portion of stray light enters the emission channel. Crisp, high-contrast TIRF images have been reported for p-TIRF [1-3, 8].

**p-TIRF stray light mitigation.** If your application permits to enclose your specimen into a closed flow cell, as shown in Fig. 4, p-TIRF geometry is the best choice. Typically, in the case of p-TIRF, the intensity of stray light is negligibly small, provided that high-quality optical materials were selected to minimize the scatter and autofluorescence inside the prism and TIRF slide. Additionally, for demanding TIRF experiments, the surface quality of the TIRF prism and TIRF slide should be high enough to ensure minimal scatter. At the emission channel, the selection of an efficient emission filter, which blocks the stray light, will reduce the background.

A limitation of p-TIRF is related to the fact that experiments with live cells in open dishes is difficult to perform with p-TIRF on an inverted microscope. An upright microscope with a prism-down geometry and water-immersion objective is a geometry well-suited for live cell studies. Visit [www.tirf-labs.com](http://www.tirf-labs.com) for additional information.

**Lightguide-based TIRF geometry.** Lg-TIRF is a sensible alternative to o-TIRF. It is a flexible geometry available with open perfusion chambers on inverted microscopes, and closed flow cells [11]. Lg-TIRF can be used with dry, water- and oil-immersion objectives, with UV and visible excitation. Novel method of Shallow Angle Fluorescence Microscopy, in which objects a few microns away from the surface are illuminated, can be implemented at the same platform [12]. Lg-TIRF principles are close to that of p-TIRF: the excitation light is naturally independent from the emission channel. In Lg-TIRF, the excitation light undergoes multiple reflections from the top and the bottom surface of glass or silica coverslip which serves as the TIRF lightguide. The excitation light escapes from the opposite end of the coverslip. Unlike in o-TIRF, it does not enter the emission channel. The optical scheme of Lg-TIRF allows for using Petri dish chambers for experiments with live cells. Potential source of scatter and autofluorescence is the coverslip material. Additionally, the top, bottom, and end surfaces of the coverslip play role in Lg-TIRF. For demanding TIRF experiments these surfaces should be of high quality.





**Ig-TIRF stray light mitigation.** To minimize the intensity of stray light and to maximize the signal-to-background ratio in Ig-TIRF, it is recommended to use coverslips made from high-quality glass or silica with low autofluorescence. The use of coverslips polished to the highest optical surface quality reduces scatter from the surfaces. However, regular coverslips work well for many Ig-TIRF applications. Higher quality coverslips as well as superior emission filters are necessary only for demanding TIRF applications.

Overall, it appears that it was not a good idea to use the emission path for delivering TIRF excitation in o-TIRF. Indeed, p-TIRF and Ig-TIRF geometries provide superior signal-to-background ratios, because the excitation light is naturally separated from the emission channel. Ig-TIRF performs all the functions that o-TIRF does. Additionally, it can be used with UV excitation. If your applications involve TIRFing of living cells in open chambers, Ig-TIRF is the optimal system for this purpose. If you are TIRFing single molecules with bright fluorophores, an inexpensive objective 1.20 NA will be sufficient. Use of expensive >1.4 NA lenses will increase the amount of light collected by the objective, and the Ig-TIRF geometry give you the freedom to select an objective, which is better fitted for your specific task. Visit [www.tirf-labs.com](http://www.tirf-labs.com) for additional information.

**CCD camera.** The TIRF signal is inherently low. Therefore, a low light CCD camera is a necessary module for TIRFing. Electron-multiplying CCD cameras, for example an iXon EMCCD by Andor Technology, is a rational choice. Scientific-grade low light CCD and sCMOS cameras, and modest performance cameras available from Thorlabs and other vendors are sufficiently sensitive for many TIRF applications, including single molecule detection [13].


**UV-absorbing-blue-light-emitting fluorophores.** p-TIRF and Ig-TIRF systems can be used with UV excitation, a feature not available in o-TIRF. There is a new class of fluorescent probes with large Stokes shift that are excited by UV light and emit fluorescence in the range 360-450 nm. These fluorophores, including modified nucleic and amino acids possess great potential for life science and material science applications. p-TIRF and Ig-TIRF systems are well suited for TIRFing such fluorophores.

**Selecting microscope objective.** Unlike in rigid o-TIRF scheme, p-TIRF and Ig-TIRF systems can be used with high and low NA objectives, dry, water- and oil-immersion lenses. This freedom causes good problem of choice. What is the best objective for your application? With p-TIRF and Ig-TIRF it entirely depends on your application. Supersensitive TIRF experiments, such as single molecule detection, require the collection of a significant portion of emitted fluorescence. The efficiency of collecting light by a microscope objective, **G**, is proportional to the numerical aperture of the objective, **NA**, to the power of 4, and the magnification, **X**, to the power of negative 2:

$$G = kS (NA)^4 (X)^{-2} \quad (4)$$

where **k**- is a coefficient, constant for the given microscope, **S**- coefficient of the CCD camera, which is typically proportional to the area of single pixel and the quantum yield of the camera. The table below enumerates the relative efficiency of objectives with different numerical apertures and magnifications. For example, an 1.20 NA X60 objective, which is available as a cost-efficient water-immersion objective, provides 42% efficiency relatively to the costly 1.49NA X60 lens and even exceeds the 1.49NAX100 objective.

The proximity of the glass surface results in anisotropy of the fluorescence emitted in TIRF mode by a fluorophore located near the surface [14, 15]. The emission patterns of fluorophores that are located closer than 50 nm to an optical interface (e.g. glass-water) are shifted toward the coverslip, which is an additional rational for using high NA objectives. However, for most single molecule detection experiments the NA=1.20 turns to be sufficient.



Compare TIRF Geometries <-> Stray Light Causes Errors in Objective-type TIRF



ev  
W  
CO  
1.  
rec  
2.  
3.  
4.  
5.  
6.  
7.  
8.  
9.  
10.  
11.  
12.  
13.  
14.  
15.



## Compare TIRF Microscopy Geometries Stray Light Errors in Objective-type TIRF

**Table.** Relative light-collecting efficiency of microscope objectives with different values of numerical aperture, NA, and magnification, X.

Numerical Aperture, NA / Magnification, X	NA 1.49 X60	NA 1.49 X100	NA 1.20 X60	NA 0.80 X40	NA 1.20 X63	NA 1.20 X100
Relative efficiency of collecting light	100%	37%	42%	19%	18%	16%

In summary, the problems of stray light, poor signal-to-background ratio, and compromised quality of the evanescent wave are inherent only to o-TIRF geometry. Significant interferences of stray light with TIRF effect - undesirable excitation of the bulk of the specimen have been documented in the literature only for the case of o-TIRF geometry. [1-6]. These interferences are caused by the fact that the excitation path uses the emission channel to deliver the excitation light to TIRF interface. On the other hand, there is substantial body of evidence, which suggests that p-TIRF and Ig-TIRF geometries demonstrate excellent quality of TIRF effect, minimum amount of stray light, crisp and high-contrast TIRF images. Both p-TIRF and Ig-TIRF geometries employ optical schemes in that the excitation light is naturally independent from the emission channel. Published articles show that p-TIRF and Ig-TIRF geometries are well-suited for single molecule detection, living cell membrane studies, and other areas of lifesciences. For supersensitive TIRF experiments it is rational to use p-TIRF and Ig-TIRF with microscope objectives of high or moderate NA numbers (>1.0 and more) to ensure that the objective collects sufficient portion of the emitted fluorescence. For additional information please visit [www.tirf-labs.com](http://www.tirf-labs.com) or contact TIRF Labs: [info@tirf-labs.com](mailto:info@tirf-labs.com).


### Literature cited:

1. Lamichhane R, Liu JJ, Pljevaljcic G, White KL, van der Schans E, Katritch V, Stevens RC, Wüthrich K, Millar DP. Single-molecule view of basal activity and activation mechanisms of the G protein-coupled receptor  $\beta$ 2AR. PNAS. 2015 Nov 17;112(46):14254-9.
2. Lamichhane R, Berezhna SY, Gill JP, Van der Schans E, Millar DP. Dynamics of site switching in DNA polymerase. J Am Chem Soc. 2013 Mar 27;135(12):4735-42.
3. Ambrose W, Goodwin P, Nolan J. Single-molecule detection with TIRF: comparing signal-to background in different geometries. Cytometry 1999, 36(3), 224.
4. Brunstein M, Teremetz M, Hérault K, Tourain C, Oheim M. Eliminating unwanted far-field excitation in objective-type TIRF. Part I. Biophys J. 2014; 106(5): 1020.
5. Brunstein M, Hérault K, Oheim M. Eliminating unwanted far-field excitation in objective-type TIRF. Part II. Biophys J. 2014; 106(5): 1044.
6. Conibear P, Bagshaw C. A comparison of optical geometries for combined flash photolysis and TIRF microscopy. J Microsc, 2000, 200(3): 218-29.
7. Mattheyses A, Axelrod D. Direct measurement of the evanescent field profile produced by objective-based TIRF. J Biomed Opt, 2006, 11: 014006A.
8. Schwarz J, König I, Anderson K. Characterizing system performance in total internal reflection fluorescence microscopy. Meth Mol Biol, 2011, 769: 373-86.
9. <https://www.semrock.com/flatness-of-dichroic-beamsplitters-affects-focus-and-image-quality.aspx>
10. <https://www.chroma.com/knowledge-resources>
11. Asanov A, Zepeda A, Vaca L. A novel form of Total Internal Reflection Fluorescence Microscopy (LG-TIRFM) reveals different and independent lipid raft domains in living cells. Biochim Biophys Acta. 2010 Feb;1801(2):147-55. See also references at [www.tirf-labs.com/applications.html](http://www.tirf-labs.com/applications.html)
12. Asanov A, Zepeda A, Vaca L. A platform for combined DNA and protein microarrays based on total internal reflection fluorescence. Sensors (Basel). 2012;12(2):1800-15.
13. Protasenko V, Hull KL, Kuno M. Demonstration of a Low-Cost, Single-Molecule Capable, Multimode Optical Microscope. Chem. Educator 2005, 10, 269282.
14. Ruckstuhl T, Enderlein J, Jung S, and Seeger S, "Forbidden Light Detection from Single Molecules," Anal. Chem. 72(9), 21172123 (2000).
15. Enderlein J, Ruckstuhl T, "The efficiency of surface-plasmon coupled emission for sensitive fluorescence detection," Opt. Express 13(22), 88558865 (2005).
16. Steyer JA, Almers W. A real-time view of life within 100 nm of the plasma membrane. Nat Rev Mol Cell Biol. 2001, 2(4), 268.





### Single ion Channel Single Molecule Detection



### iDiagnostics

cellphone based molecular diagnostics



We extended TIRF into 3<sup>rd</sup> dimension and invented iDiagnostics  
Now you can hold a hospital laboratory in the palm of your hand

### Turnkey Single Molecule Detection TIRF Microscopy System



#### Modular TIRFM systems include:


- Fluorescence microscope
- Ig-, p-, or/and o-TIRF microscopy flow systems
- Low light EM CCD camera
- Multi-color computer-controlled illuminator
- Digital fluidics SmartFlow
- Optional temperature and electric field control
- Software for instrument control and data analysis

### Prism- and Lightguide-based TIRF Microscopy Accessories

- Single molecule detection, cell membrane studies
- Superior signal-to-background ratio
- Minimal stray light, crisp, high-contrast TIRF images
- Work with dry, water-, and oil-immersion objectives
- Use UV or visible excitation light 190-900 nm
- Use Petri-dish, open perfusion, or closed flow chamber
- Nested design - fits inside 96-well plate, K-frame.  
4-inch round, or manual XY stage
- Optional temperature and electric field control



### TIRF Accessories for Fluorometers



- **TIRF Accessory TA-1004** transforms a spectrofluorometer into a super-sensitive TIRF biosensor instrument
- Optional electrochemical, DEP and temperature control
- **SmartFlow** Fluidic System allows to run unattended TIRF experiments, measure sensograms to derive  $k_{on}$  and  $k_{off}$
- Microfluidic system allows for handling nanoliter volumes